



93383

Mercuric Iodide Red extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 7774-29-0

Molecular Formula : HgI₂

Molecular Weight : 454.40

Specifications

ExiPlus	meets compendial specs of USP, Ph.Eur
Appearance (Colour)	Red
Appearance (Form)	Crystalline powder
Assay (T)	min. 99%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.0005%
Sodium (Na)	max. 0.005%
Potassium (K)	max. 0.005%
Mercurous Mercury (as Hg)	max. 0.1%
Soluble Mercury Salts (as Hg)	max. 0.05%
Solubility in KI solution	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1638 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 Gms	28521000 (GST 18%)
25 Gms	28521000 (GST 18%)
250 Gms	28521000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com