



93663

ProxiB Lambda PstI Digest Marker : 0.25-11.5kb

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
Absence of nucleases	No degradation in the marker after incubation at 37°C for 24 hours
Agarose gel electrophoresis	Well defined bands with clear intensity
Concentration	~100ng/μl

Other Information

Description

Lambda PstI Digest Marker is supplied as a ready-to-load preparation that contains 1mcg /10mcl . PstII digestion of lambda DNA generates 29 fragments (in Kb) 11.5,5.0,4.7,4.5, 2.8,2.5, 2.4, 2.4, 1.9,1.7,1.1, and 1.0 are visible. If a UV transilluminator with low sensitivity is used for detecting the DNA marker, visualizing the smallest fragments (0.22 and 0.11 kb) might require a 20 mcl load.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
50 mcg	38229090 (GST 12%)

Available Packages

50 mcg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com