



93665

Phthalic Anhydride ACS, 99%

Part A

CAS : 85-44-9

Molecular Formula : C₈H₄O₃

Molecular Weight : 148.12

Specifications

Appearance (Colour)	White
Appearance (Form)	Flakes /Solid
Solubility (Turbidity) 5% solution in Pyridine	Clear
Solubility (Colour) 5% solution in Pyridine	Colouless
Assay (NT)	min.99%
Melting Point	131° - 134°C
Residue after ignition	max.0.01%
Chloride (Cl)	max.0.002%
Sulphate (SO ₄)	max.0.003%
Heavy Metals (Pb)	max.0.0005%
Iron (Fe)	max.0.0005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2214 , IMCO Class No.:8 , Packing Group:III
HSN Code	
5 Kg	29173500 (GST 18%)
500 Gms	29173500 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com