



93762

## L-Lactic Acid pure, 88%

Part A

CAS : 79-33-4

Molecular Formula : C<sub>3</sub>H<sub>6</sub>O<sub>3</sub>

Molecular Weight : 90.08

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (T)	min. 88%
Chloride (Cl)	max. 0.01%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.001%
Sulphated Ash	max. 0.03%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 ml	29181110 (GST 18%)
500 ml	29181110 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

### Available Packages

500 ml

2500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com