



93940

Ninhydrin extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 485-47-2

Molecular Formula : C₉H₆O₄

Molecular Weight : 178.14

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Colour)	Pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 1% aq. solution	Clear
Solubility (Colour) 1% aq. solution	Pale yellow
Assay (UV)	min. 99%
Melting Point (Decomposition)	250 - 255°C
Sulphated Ash	max. 0.1%
Sensitivity to amino acids	Passes test
Suitability as Reagent for Amino acids	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29142990 (GST 18%)
100 Gms	29142990 (GST 18%)
25 Gms	29142990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

10 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com