



94408

## Sodium Thiocyanate pure, 98%

Part A

CAS : 540-72-7

Molecular Formula : NaSCN

Molecular Weight : 81.07

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (T)	min. 98%
Chloride (Cl)	max. 0.02%
Sulphate (SO4)	max. 0.01%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	28429090 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com