



94453

1,3-Dichlorobenzene extrapure, 99%

Part B

CAS : 541-73-1

Molecular Formula : C₆H₄Cl₂

Molecular Weight : 147.00

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.287-1.289
Refractive Index (20°C)	1.545-1.547
Boiling Range	171-174°C
Non Volatile Matter	max. 0.05%
Acidity (HCl)	max. 0.005%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3082 , IMCO Class No.:9 , Packing Group:III
HSN Code	
1000 ml	29039990 (GST 18%)
250 ml	29039990 (GST 18%)
2500 ml	29039990 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
250 ml	Glass Bottle

Available Packages

250 ml

1000 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com