



94601

Sodium Metasilicate Pentahydrate (Sodium Silicate Hydrated) extrapure

Part A

CAS : 10213-79-3

Molecular Formula : Na₂O₃Si₅H₂O

Molecular Weight : 212.14

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crystalline powder
Silicate (SiO ₂) content	min. 29%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3253 , IMCO Class No.:8 , Packing Group:III
HSN Code	
1 Kg	28391100 (GST 18%)
5 Kg	28391100 (GST 18%)
500 Gms	28391100 (GST 18%)
Type of Packing	
1 Kg	Plastic Bottle
5 Kg	Plastic Jar
500 Gms	Plastic Bottle

Available Packages

500 Gms
1 Kg
5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com