



94777

Propylene Glycol Monomethyl Ether 2-Acetate (MPA) extrapure, 98%

Part A

CAS : 108-65-6

Molecular Formula : C₆H₁₂O₃

Molecular Weight : 132.16

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Almost colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 25
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	0.968 - 0.970
Refractive Index (20°C)	1.401 - 1.402
Boiling Point	145 - 146°C
Acidity (CH ₃ COOH)	max. 0.05%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3271 , IMCO Class No.:3 , Packing Group:III
HSN Code	
500 ml	29091900 (GST 18%)
2500 ml	29091900 (GST 18%)

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com