



94912

1,1-Dianthrimide extrapure AR, 98%

Part B

CAS : 82-22-4

Molecular Formula : C₂₈H₁₅NO₄

Molecular Weight : 429.42

Specifications

Appearance (Colour)	Dark red to brown to orange
Appearance (Form)	Powder to crystals
Assay (HPLC)	min. 98%

Other Information

Application	Reagent for the determination of Boron.
-------------	---

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29223900 (GST 18%)
5 Gms	29223900 (GST 18%)
Type of Packing	
1 Gms	Plastic Bottle

Available Packages

1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com