



94977

Crucible Porcelain, 40ml, Tall with Lid

Part L

Specifications		
Description	White Porcelain	
Crucible	Top OD 50mm	
Outer Height	40mm	
Wall Thickness	1.8mm	
Weight	23.0gm	
Max use temperature	1150°C	
Flextural Strength @ 20°C	Min. 80 Mpa (815 Kg/cm2)	
Lid Dia	53mm	

General Information		
Storage	25 to 40°C (Room Temperature)	
IMDG Identification	Not Regulated for Transport (Non-Haz)	
HSN Code		
1 Nos	69091100 (GST 18%)	

Available Packages	
1 Nos	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com