



95056

## Potassium Sorbate pure, 99%

Part B

CAS : 24634-61-5

Molecular Formula : C<sub>6</sub>H<sub>7</sub>KO<sub>2</sub>

Molecular Weight : 150.22

### Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Powder or granules
Solubility (Turbidity) 5% aq. solution	Clear
Assay (NT)	min. 99%
Loss on drying	max. 1%
Acidity or alkalinity	Passes test
Chloride (Cl)	max. 0.005%
Iron (Fe)	max. 0.005%
Heavy Metals (Pb)	max. 0.005%
Calcium (Ca)	max. 0.005%
Copper (Cu)	max. 0.005%
Cadmium (Cd)	max. 0.005%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	29161990 (GST 18%)
500 Gms	29161990 (GST 18%)
Type of Packing	
1 Kg	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

500 Gms

1 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com