



95205

Triethylamine extrapure AR, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 121-44-8

Molecular Formula : C₆H₁₅N

Molecular Weight : 101.19

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Miscibility in Ethanol (96%)	Clear, colourless
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.726-0.728
Refractive Index (20°C)	1.400-1.402
Boiling Range	88-90°C
Non Volatile Matter	max. 0.003%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0001%
Copper (Cu)	max. 0.00001%
Water (KF)	max. 0.1%
Absorbance (A) of 6% in dilute H ₂ SO ₄ @250 nm	max. 0.04

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1296 , IMCO Class No.:3(8) , Packing Group:II
HSN Code	
1000 ml	29212990 (GST 18%)
2500 ml	29212990 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
2500 ml	Glass Bottle

Available Packages

1000 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
