



95323

## Buffer capsules pH 6 (10 capsules)

Part L

### Specifications

pH

6.0 +/- 0.05  
(Dissolve contents of 1 capsule  
& dilute upto 100 ml distilled water to  
prepare solution of pH 6 at 20°C )

### General Information

Storage

25 to 40°C (Room Temperature)

Shelf Life

60 Months

IMDG Identification

Not Regulated for Transport (Non-Haz)

HSN Code

10 caps

38229090 (GST 12%)

### Available Packages

10 caps

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com