



95365

## D-Proline extrapure CHR, 99%

Part A

CAS : 344-25-2

Molecular Formula : C<sub>5</sub>H<sub>9</sub>NO<sub>2</sub>

Molecular Weight : 115.13

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless to faint yellow
Assay	min. 99%
Specific rotation [?]	+84° to +85° (C=5 , Water)
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.1%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Homogeneity by TLC	Passes test

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29224990 (GST 18%)
5 Gms	29224990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

### Available Packages

5 Gms  
25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com