



95407

## Ammonium Bismuth Citrate extrapure AR, 43-46% Bi

Part A

CAS : 31886-41-6

Molecular Formula : C<sub>24</sub>H<sub>20</sub>Bi<sub>4</sub>O<sub>28.6</sub>NH<sub>3</sub>.10H<sub>2</sub>O

Molecular Weight : ~1875.00

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 3% aq. solution	Clear
Solubility (Colour) 3% aq. solution	Colourless
Assay	43 - 46% (Bismuth content)
pH (3% aq. solution)	5.5 - 7.5

### Other Information

Notes	Bismuth content.....43-46%
-------	----------------------------

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29181590 (GST 18%)
500 Gms	29181590 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com