



95423

## 1,2,3,4-Di-O-Isopropylidene L-Arabinopyranoside extrapure

Part B

Molecular Formula : C<sub>11</sub>H<sub>18</sub>O<sub>5</sub>

Molecular Weight : 230.16

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
10 Gms	29420090 (GST 18%)

### Available Packages

10 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com