



95425

**SWCNT Type 7 Carbon Nanotubes Single-Walled  
in aq. media pure, Length 0.5-4µm, OD 1-2nm,  
99%**

Part C

CAS : 308068-56-6

**Specifications**

Appearance (Clarity)	Hazy
Appearance (Form)	Suspension
Assay	min. 90%
OD	1 - 2 nm
Length	0.5 - 4 µm
Concentration	1mg/100 ml

**General Information**

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
20 ml	28030090 (GST 18%)

**Available Packages**

20 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com