



95536

1,4-Dioxane extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 123-91-1

Molecular Formula : C4H8O2

Molecular Weight : 88.11

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	1.032-1.034
Refractive Index (20°C)	1.421-1.423
Boiling Range	100-102°C
Freezing Point	min. 11°C
Non Volatile Matter	max. 0.001%
Acidity (CH3COOH)	max. 0.01%
Iron (Fe)	max. 0.00005%
Heavy Metals (Pb)	max. 0.0001%
Copper (Cu)	max. 0.00001%
Acetal (C6H14O2)	max. 0.3%
Peroxide (H2O2)	max. 0.005%
Carbonyl compounds (HCHO)	max. 0.01%
Water (KF)	max. 0.05%
Stabilizer	2,6 Di-tertubutylp- resol(2.5ppm)

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1165 , IMCO Class No.:3 , Packing Group:II
HSN Code	
2500 ml	29329990 (GST 18%)
500 ml	29329990 (GST 18%)
Type of Packing	

2500 ml	Glass Bottle
---------	--------------

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com