



95956

## Naringenin pure, 95%

Part B

CAS : 67604-48-2

Molecular Formula : C<sub>15</sub>H<sub>12</sub>O<sub>5</sub>

Molecular Weight : 272.26

### Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Powder
Assay (HPLC)	min.95%
Melting Point	250 - 255°C

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29143990 (GST 18%)
25 Gms	29143990 (GST 18%)
5 Gms	29143990 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

1 Gms

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com