



96045

Silver Nitrate ExiPlus, Multi-Compendial, 99.9%

Part A

CAS : 7761-88-8

Molecular Formula : AgNO₃

Molecular Weight : 169.87

Specifications

ExiPlus	meets compendial specs of USP, BP
Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (T)	min. 99.9%
Chloride (Cl)	max. 0.0005%
Free acid	Passes test
Substances not precipitated by HCl	max. 0.01%
Sulphate (SO ₄)	max. 0.0025%
Iron (Fe)	max. 0.0002%
Copper (Cu)	max. 0.0002%
Lead (Pb)	max. 0.001%
Potassium (K)	max. 0.01%
Sodium (Na)	max. 0.002%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1493 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
100 Gms	28432100 (GST 18%)
25 Gms	28432100 (GST 18%)
500 Gms	28432100 (GST 18%)

Available Packages

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com