



96240

Xylene Sulphur Free extrapure, 99%

Part A

CAS : 1330-20-7

Molecular Formula : C6H4(CH3)2

Molecular Weight : 106.17

Specifications

Appearance (Clarity)	Clear
Appearance (Form)	Liquid
Total Assay (GC, C8H10 isomers + Ethylbenzene)	min. 99%
Density (g/ml) @ 20°C	0.855 - 0.865
Refractive Index (20°C)	1.497-1.499
Boiling Range	137-142°C
Non Volatile Matter	max. 0.1%
Sulphur Compound (CS2)	max. 0.0003%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1307 , IMCO Class No.:3 , Packing Group:III
HSN Code	
25 Ltr	29024400 (GST 18%)
2500 ml	29024400 (GST 18%)
500 ml	29024400 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml
2500 ml
25 Ltr

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com