



96249

N-a-Benzoyl-DL-Arginine-4-Nitroanilide Hydrochloride (BANI, L-BAPNA) extrapure, 98%

Part B

CAS : 911-77-3

Molecular Formula : C₁₉H₂₂N₆O₄.HCl

Molecular Weight : 434.89

Specifications

Appearance (Colour)	Yellow
Appearance (Form)	Powder
Assay (T)	min. 98%

Other Information

Application	A chromogenic trypsin substrate.
-------------	----------------------------------

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29420090 (GST 18%)
100 Mg	29420090 (GST 18%)
250 Mg	29420090 (GST 18%)

Available Packages

100 Mg

250 Mg

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com