



96370

## Water, nuclease-free

Part E

CAS : 7732-18-5

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Colourless
Appearance (Clarity)	Clear
Absence of nucleases	Incubation of nucleic acids in Water, Nuclease free should show no degradation in comparison to the control as visualized after agarose gel electrophoresis.
Contamination	Sterility check kept overnight for bacterial growth

### Other Information

Description	The Water, nuclease-free, is deionized and 0.22 µm membrane-filtered. It is suitable for use in all molecular biology applications.
-------------	---

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)
5 x1 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

### Available Packages

5 x1 ml  
100 ml  
500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com