



96704

Dish Evaporating Porcelain, 100ml, OD 75mm, RB, Deep with Spout

Part L

Specifications

Appearance (Colour)	White to Off white
Appearance (Form)	Evaporating Dish
Material of Construction (MOC)	Porcelain
Evaporating Dish Top Outer Dia	75mm ± 3mm
Evaporating Dish Height	36.5mm ± 5mm
Evaporating Dish Weight	78gm ± 5gm
Max Usable Temperature	1150°C ± 50°C
Flextural Strength @20°C	Min. 80 Mpa (815 Kg/cm ²)

Other Information

Application

A high quality evaporating basin for critical laboratory use.

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
6 Nos	69091100 (GST 18%)

Available Packages

6 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com