



Sisco Research
Laboratories
Pvt. Ltd.

Product Datasheet
and Specification

96712

Chloroform extrapure, 99%

Part A

CAS : 67-66-3

Molecular Formula : CHCl3

Molecular Weight : 119.38

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC) (excluding stabilizer)	min. 99%
Density (g/ml) @ 20°C	1.474-1.480
Refractive Index (20°C)	1.444-1.445
Boiling Range	60-62°C
Non Volatile Matter	max. 0.002%
Water (KF)	max. 0.1%
Stabilizer (Ethanol)	1-2% ethanol

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1888 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
25 Ltr	29031300 (GST 18%)
2500 ml	29031300 (GST 18%)
500 ml	29031300 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com