



96763

Di-tert-Butyldicarbonate (BOC Anhydride, DiBOC) extrapure, 98.5 -101.5%

Part B

CAS : 24424-99-5

Molecular Formula : C10H18O5

Molecular Weight : 218.25

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Semi solid / Liquid
Assay (Morpholine Titration by HCl)	98.5 - 101.5 %
Melting Point	22 - 24°C

Other Information

Application

Reagent for the introduction of the Boc protection group into amino acid chains during peptide synthesis.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2930 , IMCO Class No.:6.1(4.1) , Packing Group:II
HSN Code	
100 Gms	29420090 (GST 18%)
25 Gms	29420090 (GST 18%)
5 Kg	29420090 (GST 18%)
500 Gms	29420090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages	
25 Gms	
100 Gms	
500 Gms	
5 Kg	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com