



97039

3-Trimethylsilylpropionic Acid (TMSPA) extrapure, 95%

Part B

CAS : 5683-30-7

Molecular Formula : C₆H₁₄O₂Si

Molecular Weight : 146.26

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 95%
Density (g/ml) @ 20°C	0.920-0.922
Refractive Index (20°C)	1.426-1.428

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	24 Months
HSN Code	
1 ml	29420090 (GST 18%)
5 ml	29420090 (GST 18%)

Available Packages

1 ml

5 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com