



97328

Fumaric Acid extrapure AR, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 110-17-8

Molecular Formula : C₄H₄O₄

Molecular Weight : 116.07

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in 10% NaOH	Clear
Solubility (Colour) 10% solution in 10% NaOH	Colourless
Assay (NaOH titration)	min. 99%
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.1%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	29171950 (GST 18%)
500 Gms	29171950 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com