



97628

Antimony Trichloride extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 10025-91-9

Molecular Formula : SbCl₃

Molecular Weight : 228.11

Specifications

ExiPlus	Meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Assay	min. 99%
Hydrochloric acid insoluble matter	max. 0.005%
Chloroform insoluble matter	max. 0.005%
Sulphate (SO ₄)	max. 0.05%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.003%
Copper (Cu)	max. 0.0005%
Arsenic (As)	max. 0.0005%
Sodium (Na)	max. 0.02%
Potassium (K)	max. 0.01%
Calcium (Ca)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1733 , IMCO Class No.:8 , Packing Group:II
HSN Code	
100 Gms	28273990 (GST 18%)
500 Gms	28273990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
