



97799

MES12: ICP 68 Element Multi Standard 2

Part B

Specifications	
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to Pale yellow
Appearance (Form)	Liquid
Silver (Ag)	100.0 ± 0.5 ?g/mL
Germanium (Ge)	100.0 ± 0.5 ?g/mL
Hafnium (Hf)	100.0 ± 0.5 ?g/mL
Molybdenum (Mo)	100.0 ± 0.5 ?g/mL
Niobium (Nb)	100.0 ± 0.5 ?g/mL
Antimony (Sb)	100.0 ± 0.5 ?g/mL
Silicon (Si)	100.0 ± 0.5 ?g/mL
Tin (Sn)	100.0 ± 0.5 ?g/mL
Tantalum (Ta)	100.0 ± 0.5 ?g/mL
Titanium (Ti)	100.0 ± 0.5 ?g/mL
Tungsten (W)	100.0 ± 0.5 ?g/mL
Zirconium (Zr)	100.0 ± 0.5 ?g/mL

General Information	
Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
HSN Code	
100 ml	38229090 (GST 12%)

Available Packages	
100 ml	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com