



97979

**Sulphanilic Acid (High Purity) extrapure AR,  
ACS, ExiPlus, Multi-Compendial, 99.5%**

Part A

CAS : 121-57-3

Molecular Formula : C<sub>6</sub>H<sub>7</sub>SO<sub>3</sub>N

Molecular Weight : 173.19

**Specifications**

ExiPlus	Meets compendial specs of BP
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in 5% Na <sub>2</sub> CO <sub>3</sub>	Clear
Assay (NaOH titration)	min. 99.5%
Loss on drying	max. 0.5%
Sulphated Ash	max. 0.01%
Chloride (Cl)	max. 0.002%
Sulphate (SO <sub>4</sub> )	max. 0.005%
Iron (Fe)	max. 0.005%
Heavy Metals (Pb)	max. 0.001%
Copper (Cu)	max. 0.001%
Nitrite (NO <sub>2</sub> )	max. 0.00005%
Insoluble in Sodium carbonate solution	max. 0.02%
Sensitivity for nitrite (NO <sub>2</sub> )	Passes test

**General Information**

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29214234 (GST 18%)
500 Gms	29214234 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

## Available Packages

100 Gms

500 Gms

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---