



98017

Stannic Chloride 1M in Toluene

Part A

CAS : 7646-78-8

Molecular Formula : SnCl_4

Molecular Weight:260.52

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow to yellow fuming
Appearance (Form)	Liquid
Concentration	0.95-1.10M

General Information

Storage	8 to 25°C (Cool & Dry Area)
HSN Code	
100 ml	38229090 (GST 12%)

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com