



98063

## Antimony Trioxide extrapure AR, 99%

Part A

CAS : 1309-64-4

Molecular Formula : Sb<sub>2</sub>O<sub>3</sub>

Molecular Weight : 291.50

### Specifications

Appearance (Colour)	White
Appearance (Form)	Amorphous powder
Assay (Iodometric)	min. 99%
Chloride (Cl)	max. 0.05%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.003%
Heavy Metals (Pb)	max. 0.002%
Copper (Cu)	max. 0.001%
Arsenic (As)	max. 0.005%
Calcium (Ca)	max. 0.005%
Potassium (K)	max. 0.002%
Sodium (Na)	max. 0.002%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28258000 (GST 18%)
500 Gms	28258000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com