



98087

Glutamic Oxaloacetic Transaminase (GOT) ex. Porcine Heart, 180U/mg Solids

Part B

CAS : 9000-97-9

Specifications

Appearance (Colour)	White to Yellow to Pale green
Appearance (Form)	Crystalline powder
Activity	min. 180 U/mg solids
Unit Definition	One unit will convert 1.0 micromole of alpha-ketoglutarate to L-glutamate per min at pH 7.5 at 37° C in the presence of L-aspartic acid.
Enzyme solution for activity determination	Prepared in cold distilled water.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 K.Units	35079099 (GST 18%)

Available Packages

1 K.Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com