



98465

3,3-Diaminobenzidine Tetrahydrochloride (DAB.4HCl) Buffer Substrate Solution for Peroxidase

Part B

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale brown
Appearance (Form)	Liquid
Suitability as reagent for peroxidase	Passes test

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	12 Months
HSN Code	
100 ml	38229090 (GST 12%)
Type of Packing	
100 ml	Glass Bottle

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com