



98500

Tert-Dodecylmercaptan (mix. of isomers) extrapure, 98.5%

Part B

CAS : 25103-58-6

Molecular Formula : C₁₂H₂₆S

Molecular Weight : 202.40

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to almost colourless
Appearance (Form)	Liquid
Assay (GC) (mixture of Isomers)	min. 98.5%
Density (g/ml) @ 20°C	0.840-0.860
Refractive Index (20°C)	1.460-1.461
Boiling Point	240-242°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3082 , IMCO Class No.:9 , Packing Group:III
HSN Code	
100 ml	29309099 (GST 18%)
500 ml	29309099 (GST 18%)

Available Packages

100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com