



98585

## Graphene Platelet Nanopowder (GPN Type 3), 2-10nm, 99.5%

Part C

CAS : 7782-42-5

Molecular Formula : C

Molecular Weight : 12.01

### Specifications

Appearance (Colour)	Black grey
Appearance (Form)	Powder
Carbon content	min. 99.5%
APS	15 micron
Surface area	20 - 40 m <sup>2</sup> /g
Thickness	2.0 - 10.0 nm
Bulk density	~0.10 gm/ml

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	25049090 (GST 5%)
5 Gms	25049090 (GST 5%)

### Available Packages

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com