



98680

Bathophenanthroline Disulphonate Disodium
Salt Hydrate extrapure, 99%

Part B

CAS : 52746-49-3, 98645-86-4, 53744-

Molecular Formula : C24H14N2Na2O6S2.xH2O

Molecular Weight : 536.49 (anhydrous basis)

Specifications

Appearance (Colour)	White to Orange and Faint Beige to Beige and Pink
Appearance (Form)	Powder
Assay (HPLC) (on anhydrous basis)	min. 99%
Identification (IR Spectrum)	Passes test

Other Information

Application

A water soluble reagent for determination of Iron (Fe).

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29335990 (GST 18%)
100 Mg	29335990 (GST 18%)
250 Mg	29335990 (GST 18%)

Available Packages

250 Mg
1 Gms
5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com