



98808

Hippuryl-L-Phenylalanine, 99%

Part B

CAS : 744-59-2

Molecular Formula : C18H18N2O4

Molecular Weight : 326.35

Specifications

|                     |                    |
|---------------------|--------------------|
| Appearance (Colour) | White              |
| Appearance (Form)   | Crystalline powder |
| Assay               | min. 99 %          |
| Melting Point       | 140 -143 °C        |

General Information

|                     |                                       |
|---------------------|---------------------------------------|
| Storage             | 2 to 8°C (Refrigerate)                |
| Shelf Life          | 60 Months                             |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code            |                                       |
| 100 Mg              | 29379020 (GST 18%)                    |
| 250 Mg              | 29379020 (GST 18%)                    |

Available Packages

|        |
|--------|
| 100 Mg |
| 250 Mg |

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com