



99284

## FMOC-L-Leucine extrapure, 99%

Part B

CAS : 35661-60-0

Molecular Formula : C<sub>21</sub>H<sub>23</sub>NO<sub>4</sub>

Molecular Weight : 353.40

### Specifications

Appearance (Colour)	White to Off white
Appearance (Form)	Crystalline powder
Assay (NT)	min. 99%
Melting Point	151 - 154°C

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com