



99336

Sec-Butyl Alcohol extrapure AR, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 78-92-2

Molecular Formula : C4H10O

Molecular Weight : 74.12

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Solubility	1:1 in Ethanol is clear,colourless
Assay (GC)	min 99.5%
Density (g/ml) @ 20°C	0.806-0.808
Refractive Index (20°C)	1.397-1.398
Boiling Range	98-100°C
Non Volatile Matter	max. 0.001%
Acidity (CH3COOH)	max. 0.003%
Aldehyde & Ketones (C3H7CHO)	max. 0.2%
tert-Butyl Alcohol (GC)	max. 0.1%
Isopropanol (GC)	max. 0.2%
Water (KF)	max .0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1120 , IMCO Class No.:3 , Packing Group:III
HSN Code	
2500 ml	29051490 (GST 18%)
500 ml	29051490 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com