



99629

## Dextran Sulphate Sodium Salt 500 (DSS 500) ex. Leuconostoc Sp. for molecular biology

Part B

CAS : 9011-18-1

### Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Solubility (Turbidity) 1% aq. solution	Clear
Average Molecular weight	made from Dextran 500 KDa
Chloride (Cl)	max. 0.01%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max.0.0005%
DNase, RNase, Protease	Not detected
Not for medicinal use	

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29420090 (GST 18%)
100 Gms	29420090 (GST 18%)
25 Gms	29420090 (GST 18%)
250 Gms	29420090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

10 Gms  
25 Gms  
100 Gms  
250 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com