



99752

## S. Aureus Test Kit (Pre-Weighed Media)

Part D

### Other Information

#### Applications

Used for isolation, identification and cultivation of Staphylococcus aureus.

#### Composition

(55066) Mannitol Salt broth	96.02 gm/lt.
(38907) VogelJohnson Agar Base(U/P),(I/P)	61.02 gm/lt.
(23626) Mannitol Salt Agar,Harmonized	111.02 gm/lt.
(50222) Baird Parker Agar,Base	63.00 gm/lt.

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

### Available Packages

1 Kit

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com