



99823

## 2,4-Difluoroaniline extrapure, 99%

Part B

CAS : 367-25-9

Molecular Formula : C<sub>6</sub>H<sub>5</sub>F<sub>2</sub>N

Molecular Weight : 129.11

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow to brown to purple
Appearance (Form)	Liquid
Assay	min. 99%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2941 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29214190 (GST 18%)
25 Gms	29214190 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com