



99983

## 2-Methyl Tetrahydrofuran pure, 98%

Part A

CAS : 96-47-9

Molecular Formula : C<sub>5</sub>H<sub>10</sub>O

Molecular Weight : 86.13

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	0.854-0.856
Refractive Index (20°C)	1.405-1.406
Boiling Range	78-80°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2536 , IMCO Class No.:3 , Packing Group:II
HSN Code	
100 ml	29321990 (GST 18%)
2500 ml	29321990 (GST 18%)
500 ml	29321990 (GST 18%)
Type of Packing	
100 ml	Glass Bottle
2500 ml	Glass Bottle
500 ml	Glass Bottle

### Available Packages

100 ml

500 ml

2500 ml

25 Ltr

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com