



10674

## Ethyl Iodide extrapure, 99%

Part B

CAS : 75-03-6

Molecular Formula : C<sub>2</sub>H<sub>5</sub>I

Molecular Weight : 155.97

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Colour (APHA)	max. 30
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.920-1.930
Refractive Index (20°C)	1.512-1.514
Boiling Range	71-73°C
Stabilizer	silver wire

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	29035990 (GST 18%)
25 ml	29035990 (GST 18%)
250 ml	29035990 (GST 18%)
Type of Packing	
100 ml	Glass Bottle
250 ml	Glass Bottle

### Available Packages

25 ml

100 ml

250 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com