



14215

Polyanetholesulphonic Acid Sodium Salt (Sodium Polyanetholesulphonate, SPS) extrapure

Part B

CAS : 55963-78-5

Molecular Formula : (C₁₀H₁₁NaO₄S)_n

Specifications

Appearance (Colour)	White to yellow to beige to cream
Appearance (Form)	Powder
Solubility (Turbidity) 1% aq. solution	Clear
pH (1% aq. solution)	4.5 - 7.0
Loss on drying	max.15.0%
Absorbance (A) of 1% aq. solution in a 1cm cell v/s H ₂ O @283 -286 nm	min.75

Other Information

Applications

Used as an anticoagulant and inhibits antimicrobial systems in blood.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29309099 (GST 18%)
25 Gms	29309099 (GST 18%)
5 Gms	29309099 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com