



17083

Phenyl Isothiocyanate (PITC) pure, 98%

Part B

CAS : 103-72-0

Molecular Formula : C₇H₅NS

Molecular Weight : 135.19

Specifications

| | |
|-------------------------|-----------------|
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Brownish yellow |
| Appearance (Form) | Liquid |
| Assay (GC) | min. 98% |
| Density (g/ml) @ 20°C | 1.130 - 1.132 |
| Refractive Index (20°C) | 1.649 - 1.650 |

General Information

| | |
|---------------------|--|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:2927 , IMCO Class No.:6.1(8) , Packing Group:II |
| HSN Code | |
| 100 ml | 29309099 (GST 18%) |
| 500 ml | 29309099 (GST 18%) |
| Type of Packing | |
| 500 ml | Glass Bottle |

Available Packages

25 ml

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com