



18437

Bergamot Oil extrapure

Part B

CAS : 8007-75-8

Specifications

| | |
|-------------------------|--|
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Faint yellow to dark yellow green liquid |
| Appearance (Form) | Liquid |
| Odour | Characteristic odour |
| Refractive Index (20°C) | 1.461 - 1.468 |
| Density (g/ml) @ 20°C | 0.860 - 0.880 g/ml |

General Information

| | |
|---------------------|---|
| Storage | 8 to 25°C (Cool & Dry Area) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:1993 , IMCO Class No.:3 , Packing Group:III |
| HSN Code | |
| 1000 ml | 33012540 (GST 18%) |
| 500 ml | 33012540 (GST 18%) |

Available Packages

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com